



Rotor of Micro Centrifuge
6 tubes of 1.5ml/2ml capacity



Rotor of Micro Centrifuge
12 tubes of 0.5ml capacity



Model: T-1

Mini Centrifuge Digital

- Micro processor controlled, Micro-Centrifuge machine (Mini Centrifuge)
- Speed range from 1000 to 10000 r.p.m. in step of 1000.

Technical Specifications

- ▶ Fitted with maintenance free motor (Brushless without carbon)
- ▶ Complete body made of ABS plastic injection moulded.
- ▶ Fitted with 0-15 minute digital timer.
- ▶ Supply completed with 2 standard rotors as under :
 - (i) Angle Rotor 6 micro centrifuge of 1.5ml capacity
 - (ii) Angle rotor 12 micro centrifuge tubes of 0.5ml capacity

Size	Weight
R.P.M.	RCF (x g)
1000	56
2000	223
3000	503
4000	894

Size	Weight
R.P.M.	RCF (x g)
5000	1396
6000	2011
7000	2737
8000	3575

Size	Weight
R.P.M.	RCF (x g)
9000	4524
10000	5585

Mini Centrifuge Non Digital Max. Speed 7000 r.p.m.



Rotor of Micro Centrifuge
6 tubes of 1.5ml/2ml capacity



Rotor of Micro Centrifuge
12 tubes of 1.5ml/2ml capacity



Model: T-2

Mini Centrifuge (Non-Digital)

- Micro-Centrifuge machine (Mini Centrifuge).
- Speed fixed 7000 rpm.

Technical Specifications

- ▶ Fitted with high quality motor with carbon brushes.
- ▶ Complete body made of ABS plastic injection moulded.
- ▶ Supply completed with 2 standard rotors as under
 - (i) Angle Rotor 6 micro centrifuge of 1.5ml capacity.
 - (ii) Angle rotor 12 micro centrifuge tubes of 0.5ml capacity.