


Paper Electrophoresis Horizontal Vertical (3 Strips Capacity)



Model: PEH-1
Horizontal having 3 strips capacity.



Model: PEV-2
Vertical having 3 strips capacity

Specifications

- The buffer chamber is made out of single moulded piece to avoid leakage of buffers.
- It is made of transparent poly carbonate.
- The unit have 2 Platinum electrodes, supplied with their banana sockets defined in black and red colour.
- Horizontal paper strip holder to adopt 3 pare strips is supplied with horizontal unit.
- Vertical strip holder to adopt 3 paper strip in vertical position is supplied with vertical unit.
- Supply complete with respective protection cover, red and white cables.
- To work on 220 volts 50Hz A.C.
- Packed in corrugated box with instruction manual.