



Model: VM-20

Specifications

- Compact, rugged construction.
- Heavy metallic base fitted with rubber vacuum feets which prevents the movement of shaker during use.
- Choice of continuous & touch mode.
- Variable speed control.
- Maximum speed 3000 RPM.
- To Work on 220Volts_s + 10% 50 HZ A.C.
- Include following standard set up with variety of inserts.

Supply as above and completed with following accessories

▶ Rubber Cup Attachment for Test Tubes	1 No.
▶ Insert for Micro Tubes of 0.2ml/0.5ml	1 No.
▶ Insert for Micro Tubes of 1.5ml/2ml	1 No.
▶ Insert Plain for Conical/Volumetric Flask	1 No.

- Packed in corrugated box.