



Model:- EP-1

Electrophoresis Power Supply, Analog Fixed current

Specifications

- Analog Model with fixed current 250mA and also 5 types of fixed voltage means 50V, 100V, 150V, 200V and 250V.
- It is recommended for mini Horizontal and mini vertical Electrophoresis Apparatus.
- The voltage is selected thru a 5 step voltage regulator. Example if selector switch is at 1 position. The LED light will glow at 50V means fixed voltage of 50V is selected. Like wise other voltage can be selected.



Model:- EP-2

Electrophoresis Power Supply, Digital, variable

Specifications

- Body made of ABS plastic injection moulded and is shock proof.
- Microprocessor based, programmable, Electrophoresis power supply.
- Variable voltage 0-300V.
- Variable current 0-500mA.
- Digital timer 9 hr. 99 minutes programmable
- Digital display of voltage and current thru 2 line 22 character LCD display.
- Digital display of timer in respect to set time and run time.
- For single unit use at a time.
- To work on 220 volts 50 Hz A.C.
- Supply is complete with cord & plug, Instruction manual, Packed in corrugated box.